 University of Zurich <small>UZH</small> Institute of Laboratory Animal Sciences	Standard Operating Procedure SOP	Page 2 of 4
Date: 21.08.2019	CSF sampling via cisterna magna	LTK-TRT-42-A-EN Version: A

Method: Sampling of cerebrospinal fluid (CSF)

Min/Max amount:

The maximum volume of CSF collected is 8 µl / 20 g mouse. Serial sampling is possible every 2 months.

Storage of Material:


Collection tubes (Eppendorf) are found in the animal room.

Material:

1. Scalpel, straight edge razor,
2. Glass capillary,
3. Capillary puller,
4. Shaver,
5. Betadine iodine solution,
6. Cotton swabs,
7. Small animal stereotactic frame,
8. Vitamin A eye ointment,
9. Dissection microscope,
10. Surgical glue.

Safety:

1. General rules for working with sharp tools (scalpels, syringes, scissors) have to be followed.
2. Follow the rules of the animal house

 <p>University of Zurich UZH Institute of Laboratory Animal Sciences</p>	<p align="center">Standard Operating Procedure</p> <p align="center">SOP</p>	<p align="center">Page 3 of 4</p>
<p>Date: 21.08.2019</p>	<p align="center">CSF sampling via cisterna magna</p>	<p align="center">LTK-TRT-42-A-EN Version: A</p>

Method Description:

1. Pull a glass capillary and trim it so the inner diameter is about 0.5 mm^[JvB1].
2. Anesthetize the mice according to SOP LTK-TRT-18-D-EN Injection anesthesia.
3. Apply vitamin A ointment on the eyes.
4. Shave the neck of the mouse, clean it from cut hair and disinfect using iodine solution.
5. Fix mouse on a suitable stereotactic head holder so the head forms a 135° angle with the body.
6. Using a scalpel make a 1 cm sagittal incision of the skin, inferior to the occiput.
7. Under the dissection microscope separate the subcutaneous tissue. The dura mater of cisterna magna should be visible as on the figure 1.





Figure 1. Dura mater of the cisterna magna exposed through the skin incision on the neck. Arrow points the point to be penetrated with the capillary (adapted from: 1)

8. Dry the dura mater with cotton swab.
9. Penetrate the capillary tube into the cisterna magna through the dura mater as on the figure 1. Capillary is in the cisterna magna once the tissue resistance changes. CSF will flow into the capillary.
10. Remove the capillary and flush the CSF with a syringe into a clean tube.
11. Suture the skin with surgical glue.
12. Apply analgesia according to SOP-LTK-TRT-17-B-EN Post-surgery analgesia
13. Apply antidote (SOP-LTK-TRT-18-D-EN Injection anesthesia)
14. Move animal into wake up cage (a regular cage placed on a 37°C electrical heating mat, covered with a surgical cloth), only put fully awake animals back to the housing cage.
15. Check for postoperative complications after 1-2 h and the next day and re-apply analgesia if necessary (SOP-LTK-TRT-17-B-EN Post-surgery analgesia).

 University of Zurich <small>UZH</small> Institute of Laboratory Animal Sciences	Standard Operating Procedure SOP	Page 4 of 4
Date: 21.08.2019	CSF sampling via cisterna magna	LTK-TRT-42-A-EN Version: A

Documentation:

Lab Journal.

Problem management:

In case of serious adverse events, contact supervisor, lab head or vet.

Sample storage:

Cerebrospinal fluid samples need to be frozen immediately after collection.

Literature:

1. DeMattos RB, Bales KR, Parsadanian M, O'Dell MA, Foss EM, Paul SM, Holtzman DM; Plaque-associated disruption of CSF and plasma amyloid- β (A β) equilibrium in a mouse model of Alzheimer's disease. *Journal of Neurochemistry*. 2002, 81, 229-236
2. Liu L, Herukka SK, Minkeviciene R, van Groen T, Tanila H; Longitudinal observation on CSF Abeta42 levels in young to middle-aged amyloid precursor protein/presenilin-1 doubly transgenic mice. *Neurobiol Dis*. 2004 Dec;17(3):516-23.
3. Liu L, Duff K; A Technique for Serial Collection of Cerebrospinal Fluid from the Cisterna Magna in Mouse. *J. Vis. Exp*. 2008 (21), e960.